| | |
|----|--|
| UK | OXFORD PARTIAL KNEE SYSTEM ATTENTION OPERATING SURGEON |
| CZ | PARCIÁLNÍ KOLENNÍ SYSTÉM OXFORD UPOZORNĚNÍ PRO OPERUJÍCÍHO CHIRURGA |
| DE | OXFORD TEILKNEISYSTEM INFORMATION FÜR DEN OPERIERENDEN CHIRURGEN |
| DK | OXFORD PARTIELT KNÆSYSTEM BEMÆRKNING TIL OPERERENDE KIRURG |
| EL | ΣΥΣΤΗΜΑ ΜΕΡΙΚΗΣ ΑΡΘΡΟΠΛΑΣΤΙΚΗΣ ΓΟΝΑΤΟΣ OXFORD ΠΡΟΣ ΤΟΝ ΧΕΙΡΟΥΡΓΟ |
| ES | SISTEMA DE RODILLA PARCIAL OXFORD INFORMACIÓN PARA EL CIRUJANO |
| FI | OXFORD-OSITTAINEN POLVIJÄRJESTELMÄ LEIKKAVALLE KIRURGILLE |
| FR | SYSTÈME DE GENOU PARTIEL OXFORD À L'INTENTION DU CHIRURGIEN |
| HU | OXFORD RÉSZLEGES TÉRDPROTÉZISRENDSZER AZ OPERÁLÓ SEBÉSZ FIGYELMÉBE |
| IT | OXFORD SISTEMA DI GINOCCHIO PARZIALE ALL'ATTENZIONE DEL CHIRURGO |
| KO | OXFORD 부분식 슬관절 시스템 시술 시 주의 사항 |
| NL | OXFORD PARTIËLE-KNIEVERVANGINGSSYSTEEM INFORMATIE VOOR DE BEHANDELEND CHIRURG |
| NO | OXFORD DELKNESYSTEM INFORMASJON FOR KIRURGEN |
| PL | SYSTEM CZĘŚCIOWYCH ENDOPROTEZ STAWU KOLANOWEGO OXFORD UWAGI PRZEZNACZONE DLA OPERUJĄCEGO CHIRURGA |
| PT | SISTEMA DE JOELHO PARCIAL OXFORD À ATENÇÃO DO CIRURGIÃO |
| RU | СИСТЕМА ЧАСТИЧНОГО ЭНДОПРОТЕЗИРОВАНИЯ КОЛЕННОГО СУСТАВА OXFORD ВНИМАНИЮ ОПЕРИРУЮЩЕГО ХИРУРГА |
| SE | OXFORD PARTIELLA KNÄSYSTEM INSTRUKTIONER TILL OPERERANDE KIRURG |
| TR | OXFORD KISMİ DİZ SİSTEMİ CERRAHİN DİKKATİNE |
| ZH | OXFORD部分膝关节系统(OXFORD PARTIAL KNEE SYSTEM) 主刀外科医生需注意 |

0086

BIOMET UK LTD.
WATERTON INDUSTRIAL ESTATE
BRIDGEND CF31 3XA, UK
TEL: +44 (0) 1656 655221
FAX: +44 (0) 1656 645454

REF: 5401000201
REVISION: L
DATE: JUNE 2015
Change Control No:
ECR 227884

For additional information please see www.biomet.com

UK

OXFORD PARTIAL KNEE SYSTEM ATTENTION OPERATING SURGEON

DESCRIPTION

The Oxford Partial Knee is a medial partial knee replacement system consisting of a femoral component, a tibial component and a freely mobile meniscal bearing.

Materials

Femoral Component: CoCrMo Alloy
Tibial Component: CoCrMo Alloy
Meniscal Bearing: Ultra-High Molecular Weight Polyethylene (UHMWPE)
Marker Wire: Titanium Alloy
Marker Balls: Tantalum

INDICATIONS

The Oxford Partial Knee is intended for use in individuals with osteoarthritis or avascular necrosis limited to the medial compartment of the knee and is intended to be implanted with bone cement.

CONTRAINDICATIONS

Contraindications include:

1. Infection, sepsis, and osteomyelitis.
2. Use in the lateral compartment of the knee.
3. Rheumatoid arthritis or other forms of inflammatory joint disease.
4. Revision of a failed prosthesis, failed upper tibial osteotomy or post-traumatic arthritis after tibial plateau fracture.
5. Insufficiency of the collateral, anterior or posterior cruciate ligaments which would preclude stability of the device.
6. Full thickness damage to the weight bearing area of the lateral compartment of the knee.
7. Uncooperative patient or patient with neurologic disorders who are incapable of following directions.
8. Insufficient bone stock to provide adequate support and/or fixation to the prosthesis.
9. Metabolic disorders which may impair bone formation.
10. Osteomalacia.
11. Distant foci of infections which may spread to the implant site.
12. Rapid joint destruction, marked bone loss or bone resorption apparent on roentgenogram.
13. Vascular insufficiency, muscular atrophy, neuromuscular disease.
14. Incomplete or deficient soft tissue surrounding the knee.
15. Charcot's disease.
16. A fixed varus deformity (not passively correctable) of greater than 15 degrees.
17. A fixed flexion deformity (not passively correctable) of greater than 15 degrees.

WARNINGS

1. Improper selection, placement, positioning, alignment and fixation of the implant components may result in unusual stress conditions which may lead to subsequent reduction in the service life of the prosthetic components.
2. Improper preoperative or intraoperative implant handling or damage (scratches, dents, etc.) can lead to crevice corrosion, fretting, fatigue fracture and/or excessive wear.
3. Do not modify implants.
4. Do not reuse implants. While an implant may appear undamaged, previous stress may have created imperfections that would reduce the service life of the implant. Do not treat patients with implants that have been, even momentarily, placed in a different patient. Furthermore, re-using an implant could cause patient contamination.
5. Malalignment or soft tissue imbalance can place inordinate forces on the components which may cause excessive wear to the patellar or tibial bearing articulating surfaces. Revision surgery may be required to prevent component failure.
6. Care is to be taken to ensure complete support of all parts of the device embedded in bone cement to reduce risk of stress concentrations, which may lead to failure of the procedure. Complete preclosure cleaning and removal of bone cement debris, metallic debris, and other surgical debris at the implant site is critical to minimize wear of the implant articular surfaces. Implant fracture and loosening due to cement failure has been reported.
7. The surgeon is to be thoroughly familiar with the implants and instruments, prior to performing surgery.
8. Patients should be warned of the impact of excessive loading that can result if the patient is involved in an occupation that includes substantial walking, running, lifting, or excessive muscle loading due to weight that place extreme demands on the knee and can result in device failure or dislocation.
9. Patient smoking may result in delayed healing, non-healing and/or compromised stability in or around the placement site.
10. The surgical technique should be followed. Deviations from the surgical technique could result in early loosening/failure of the device, or other adverse events as outlined in the following section. Clinical outcome may be affected by component positioning. Proper placement of the implant should take into consideration individual patient anatomy as well as surgeon preference. The surgical technique sets forth guidelines for placement of the knee system.

PRECAUTIONS

- As with other surgical procedures, errors of technique are most likely when the method is being learned. To reduce these to a minimum, surgeons are required in the United States and strongly recommended throughout the world, to attend an Instructional Course on the Oxford Partial Knee before attempting the operation.
- Biomet joint replacement prostheses provide the surgeon with a means of reducing pain and restoring function for many patients. While these devices are generally successful in attaining these goals they cannot be expected to withstand the activity levels and loads of normal healthy bone and joint tissue.
- Accepted practices in postoperative care are important. Failure of the patient to follow postoperative care instructions involving rehabilitation can compromise the success of the procedure. The patient is to be advised of the limitations of the reconstruction and the need for protection of the implants from full load bearing until adequate fixation and healing have occurred. Excessive, unusual and/or awkward movement and/or activity, trauma, excessive weight, and obesity have been implicated with premature failure of the implant by loosening, fracture, dislocation, subluxation and/or wear. Loosening of the implants can result in increased production of wear particles, as well as accelerate damage to bone making successful revision surgery more difficult. The patient is to be made aware and warned in advance of general surgical risks, possible adverse effects as listed, and to follow the instructions of the treating physician including follow-up visits.
- Specialised instruments are designed for Biomet joint replacement systems to aid in the accurate implantation of the prosthetic components. The use of instruments or implant components from any other systems can result in inaccurate fit, sizing, excessive wear and device failure.
- Intraoperative fracture or breaking of instruments has been reported. Surgical instruments are subject to wear with normal usage. Instruments, which have experienced extensive use or excessive force, are susceptible to fracture. Surgical instruments should only be used for their intended purpose. Biomet recommends that all instruments be regularly inspected for wear and disfigurement.
- All trial, packaging, and instrument components must be removed prior to closing the surgical site. Do not implant.

POSSIBLE ADVERSE EFFECTS

A time-course distribution of the adverse events reported in the clinical investigation of the Oxford Partial Knee using a standard open surgical technique is provided in Table 1.

Table 1 – Time-Course Distribution of Adverse Events reported in the clinical trial for the Oxford Partial Knee* using a standard open surgical technique

| Adverse Events | | Frequency | | | | | | Percent of Population ¹ (n = 125) |
|--|--|-----------|------|------|------|------|------|---|
| Visit | | IO | 6 mo | 1 yr | 2 yr | 3 yr | 4 yr | >5 yr |
| Local - Operative Site | | | | | | | | |
| Effusion | | | | 1 | | | | 0.8% |
| Deep Infection | | | | | 1 | | | 0.8% |
| Degeneration of contralateral condyle | | | | | | 1 | 3 | 3.2% |
| Loose body and/or osteophyte removal | | | 1 | | 2 | | 1 | 3.2% |
| Soft tissue damage | | 2 | | | | | | 1.6% |
| <i>Dislocation</i> | | | | | 2 | | | 1.6% |
| <i>Component mal-alignment</i> | | 1 | | | | | | 0.8% |
| <i>Patella dislocation</i> | | | | | 1 | | | 0.8% |
| <i>Component loosening</i> | | | | | | 1 | 2 | 4.8% |
| Post-operative bone fracture | | 1 | | | | | | 0.8% |
| Trauma | | | 1 | | | | | 0.8% |
| Mechanical symptoms | | | | 1 | | | | 0.8% |
| <i>Instability</i> | | | | | | | 1 | 0.8% |
| <i>Persistent pain</i> | | | | | | 1 | | 0.8% |
| <i>Wear of bearing due to osteophyte</i> | | | | | | 1 | | 0.8% |
| <i>Systemic</i> | | | | | | | | |
| Development of rheumatoid arthritis | | | | 1 | | | | 0.8% |

* Phase 2 device design

IO = intraoperatively

¹All percentages for adverse events are based the number of occurrences reported in a patient population of 125 knee cases.

Those events listed in italics are considered device related events.

Boldface numbers represent revisions due to the given adverse event. One additional case was revised at 130 months post-operatively, cause unknown.

The following complications have also been reported in the clinical literature for partial and total knee arthroplasty designs and could potentially occur with the Oxford Partial Knee device.

- Major surgical risks associated with anaesthesia including brain damage, pneumonia, blood clots, heart attack, and death.
- Cardiovascular disorders including venous thrombosis, pulmonary embolism, and myocardial infarction.
- A sudden drop in blood pressure intraoperatively due to the use of bone cement.
- Damage to blood vessels, hematoma, delayed wound healing and/or infection.
- Temporary or permanent nerve damage may result in pain and numbness.

- Material sensitivity reactions. Implantation of foreign material in tissues can result in histological reactions involving various sizes of macrophages and fibroblasts. The clinical significance of this effect is uncertain, as similar changes may occur as a precursor to or during the healing process.
- Particulate wear debris and discoloration from metallic and polyethylene components of joint implants may be present in adjacent tissue or fluid. It has been reported that wear debris may initiate a cellular response resulting in osteolysis or osteolysis may be a result of loosening of the implant.
- Early or late postoperative, infection, and allergic reaction.
- Intraoperative bone perforation or fracture may occur, particularly in the presence of poor bone stock caused by osteoporosis, bone defects from previous surgery, bone resorption, or while inserting the device.
- Loosening or migration of the implants can occur due to loss of fixation, trauma, malalignment, bone resorption, excessive activity.
- Periarticular calcification or ossification, with or without impediment of joint mobility.
- Inadequate range of motion due to improper selection or positioning of components.
- Dislocation and subluxation due to inadequate fixation and improper positioning. Muscle and fibrous tissue laxity can also contribute to these conditions.
- Fatigue fracture of components can occur as a result of loss of fixation, strenuous activity, malalignment, trauma, non-union, or excessive weight.
- Fretting or crevice corrosion can occur at interfaces between components.
- Wear and/or deformation of articulating surfaces.
- Valgus-varus deformity.
- Transient peroneal palsy secondary to surgical manipulation and increased joint movement has been reported following knee arthroplasty in patients with severe flexion and valgus deformity.
- Patellar tendon rupture and ligamentous laxity.
- Persistent pain.

PATIENT SELECTION

Positive selection factors to be considered include:

- ACL and PCL functionally intact.
- Cartilage and bone erosions limited to the anterior and middle parts of the medial compartment. The posterior part of the medial compartment and the lateral compartment having cartilage of normal thickness.
- Medial collateral ligament not structurally shortened (i.e. varus deformity correctable).
- Patellofemoral joint damage limited to (or greater on) the medial facets.
- Fixed flexion deformity of less than 15 degrees.
- Flexion possible to 110 degrees under anaesthetic.
- Need to obtain pain relief and improve function.
- Ability and willingness of the patient to follow instructions, including control of weight and activity level.
- A good nutritional state of the patient.
- The patient must have reached full skeletal maturity.

CLINICAL STUDIES

A prospective multi-site clinical investigation of the Oxford Partial Knee involving 125 knee devices in 107 patients (see Tables 2 and 3) was conducted in the United States to determine the safety and effectiveness of the device when implanted using a standard open surgical technique. All clinical results and adverse events for this investigation were derived from the Oxford Partial Knee Phase 2 device, a previous version of the current Phase 3 device, that had a single femoral component size, a universal (left and right) design tibial component of few sizes, and a universal design meniscal bearing component with extended sizes.

Table 2 – Patient Demographics for the Oxford Clinical Study (Phase 2 Device)

| | All Oxford Knees Enrolled |
|-------------------------------|---|
| Total # Knees (# Patients) | 125 (107) |
| Mean Age in years (range) | 63±10.6 (29-85) |
| Sex | Males - 60 Females - 65 |
| Indications | Osteoarthritis - 114 Post-Traumatic Arthritis - 10 Avascular Necrosis - 1 |
| Side | Left - 56 Right - 69 |
| Compartment | Medial - 119 Lateral - 6 |
| Mean Height in Inches (range) | 67.0±3.9 (59-77) |
| Mean Weight in Pounds (range) | 187±38.6 (105-256) |

Table 3 – Device Accounting for the Oxford Clinical Study (Phase 2 Device) based on number of completed clinical follow-up examinations

| | 6 months | 1 year | 2 year | 3 year | 4 year | 5 year |
|---------------------------------|----------|--------|--------|--------|--------|--------|
| ¹ Theoretically Due | 125 | 125 | 125 | 113 | 102 | 84 |
| ² Deaths | 0 | 0 | 1 | 2 | 2 | 2 |
| ³ Revisions | 3 | 4 | 8 | 11 | 13 | 15 |
| ⁴ Expected | 122 | 121 | 116 | 100 | 87 | 67 |
| ⁵ Clinical Follow-Up | 100 | 110 | 80 | 83 | 69 | 51 |
| ⁶ Percent Follow-up | 82.0% | 90.9% | 69.0% | 83.0% | 79.3% | 76.1% |

¹ Based on the cut-off date when the last patient enrolled reached their 2 year post-operative anniversary

² Cumulative over time

³ Any component removed, cumulative over time

⁴ Theoretically Due - (Deaths + Revised)

⁵ Cases with complete clinical data (i.e., HSS, radiographic), obtained at the specified time point

⁶Clinical Follow-Up / Expected

Each patient was evaluated pre-operatively, and at the immediate and 6, 12, and 24 month post-operative intervals, and annually thereafter until the last patient enrolled had achieved their 24 month follow-up. All operative and post-operative complications, device related or not, were recorded for patients enrolled into the investigation (see Table 1).

Clinical results were evaluated using the Hospital for Special Surgery (HSS) knee scoring system and radiographic data. At each follow-up visit an HSS knee score and anterior/posterior and lateral radiographs were obtained. Radiographs were reviewed by the implanting. See Table 4 for clinical study results.

A patient was defined as a success if they met all of the following criteria:

1. A Good/Excellent HSS score, i.e., > 70 points.
2. No radiolucent lines > 1 mm in width surrounding > 50% of the component after 1 year in-situ.
3. No progressive radiolucencies.
4. No revision/removal of any components.

Table 4 – Oxford Clinical Study Results* (Phase 2 Device) using a standard open surgical technique

| | Preop | 1 year | 2 years | 3 years | 4 years | 5 years |
|---|--------------------|--------------------|------------------|------------------|------------------|------------------|
| Cases with complete HSS | 123 | 110 | 80 | 83 | 69 | 51 |
| Average HSS Score | 59.5 | 89.3 | 90.0 | 90.6 | 90.7 | 90.4 |
| ¹ Cases Rated as Good/Excellent HSS | 20/123 (16.3%) | 105/110 (95.5%) | 77/80 (96.3%) | 82/83 (98.8%) | 64/69 (92.8%) | 50/51 (98.0%) |
| Femoral Lucencies = 1 mm | 6/108 (5.5%) | 2/80 (2.4%) | 2/83 (2.4%) | 2/68 (2.9%) | 2/51 (2.9%) | |
| Tibial Lucencies = 1 mm | 5/108 (4.6%) | 6/80 (7.5%) | 8/83 (9.6%) | 7/68 (10.3%) | 3/51 (5.9%) | |
| Number of G/E cases with radiolucent lines > 1 mm around > 50% of component | 0 | 0 | 0 | 0 | 1 (femoral) | |
| Number of G/E cases with progressive radiolucencies | 0 | 0 | 0 | 1 (tibial) | 0 | |
| ² Revisions | 4 | 8 | 11 | 13 | 15 | |
| ³ Cumulative Survivorship | 96.75% | 93.34% | 90.73% | 88.83% | 86.82% | |
| ⁴ Successful Cases | 105 | 77 | 82 | 63 | 49 | |
| ⁵ Percent Successful | 92.5% (105/114) | 87.5% (77/88) | 87.2% (82/94) | 76.8% (63/82) | 74.2% (49/66) | |

* Based on the cut-off date when the last patient enrolled reached their 2 year post-operative anniversary

¹Hospital for Special Surgery score > 70

² Number of components removed at specified time point

³ Kaplan-Meier Life Table results

⁴ A successful case required a Good/Excellent HSS score, no revision/removal of any component, no radiolucent lines > 1 mm in width surrounding > 50% of the component, and no progressive radiolucencies.

⁵ Denominator includes cases with complete HSS and radiographic data, and revisions.

There were a total of 23 revisions reported for the Oxford study group (over a follow-up period of at least 9 years), with 8 of these occurring within the first 2 years of implantation. Of the 8 revisions reported at 2 years, 2 were for tibial bearing dislocation, 1 for patellar dislocation, 1 for infection, 1 for component malalignment, 1 for recurrent arthritis due to trauma, 1 for onset rheumatoid arthritis, and 1 for femoral loosening and fracture at the bone-cement interface. In all but 1 case the knees were revised to a total knee prosthesis. For the remaining 15 revisions reported after 2 years, 6 were due to loosening, 4 to progression of osteoarthritis in the lateral compartment, 1 to persistent pain, 1 to instability, 1 to impingement on an osteophyte and subsequent wear of the tibial bearing, 1 to impingement of an osteophyte on the femur, and 1 failed to report a reason. Revisions in this latter group occurred from 2 to 12 years post-operatively.

The survival rate for the Oxford Partial Phase 2 device study group at 2 years post-operatively is 93.38%, based on the endpoint of revision/removal of any component. Table 5 displays the Kaplan-Meier life table for survivorship through 8 years post-operatively for the Oxford study group. Survivorship rates for the study group are comparable to those rates seen in the literature for other partial knee devices and the rates seen in other studies of the Oxford Phase 2 device.

Table 5 – Survivorship for Oxford Clinical Study (Phase 2 Device)

| Interval Since Operation (years) | Number in Beginning of Interval | Number of Revisions at End of Interval | % Interval Survival ¹ | % Cumulative Survival ² | 95% Confidence Interval |
|----------------------------------|---------------------------------|--|----------------------------------|------------------------------------|-------------------------|
| 0-1 | 125 | 4 | 96.75% | 96.75% | (93.61 - 99.98) |
| 1-2 | 117 | 4 | 96.52% | 93.38% | (88.95 - 97.82) |
| 2-3 | 109 | 3 | 97.16% | 90.73% | (88.50 - 95.95) |
| 3-4 | 99 | 2 | 97.91% | 88.83% | (83.08 - 94.57) |
| 4-5 | 90 | 2 | 97.74% | 86.82% | (80.57 - 93.07) |
| 5-6 | 85 | 0 | 100% | 86.82% | (80.57 - 93.07) |
| 6-7 | 65 | 3 | 94.92% | 82.41% | (75.21 - 89.60) |
| 7-8 | 50 | 1 | 97.87% | 80.65% | (73.35 - 87.95) |

¹ Percent survival for that interval only, taken at the end of the interval.

² Percent cumulative survival taken at the end of the interval.

In addition, 2 year clinical data from 328 knee cases implanted with the current Phase 3 device, implanted using the minimally invasive surgical technique and minimally invasive surgical instruments specifically developed for the Phase 3 device, was collected from 3 European sites (2 U.K., 1 Holland). European clinical results were evaluated using the Knee Society Score (KSS) scoring system. At 2 years following surgery 5 of the 307 knees (1.6%) with available data had been revised (see Table 6).

Table 6 – Results at 2 years for Phase 2 Device using an open surgical technique and Phase 3 Device using a minimally invasive surgical technique

| Clinical Parameters | Oxford Study Phase 2 N = 125 knees | Combined European Data* Oxford Phase 3 N = 328 knees | European Site 1 Oxford Phase 3 N = 208 knees | European Site 2 Oxford Phase 3 N = 40 knees | European Site 3 Oxford Phase 3 N = 80 knees |
|--|---|--|--|---|---|
| Revision Rate ¹ | 6.8% (8/117) | 1.6% (5/307) | 2.0% (4/196) | 2.7% (1/37) | 0% (0/74) |
| Percent with a Good or Excellent Knee Score ² | N = 80 96.3%³ (77/80) | N = 271 83.0%⁴ (225/271) | N = 160 83.1%⁴ (133/160) | N = 37 86.5%⁴ (32/37) | N = 74 81.0%⁴ (60/74) |

* Combined data from European Site 1, Site 2, and Site 3.

European Site 1 = Nuffield Orthopaedic Centre (U.K.), Site 2 = Macclesfield Hospital (U.K.), and Site 3 = Groningen Hospital (Holland).

¹ Revision rate (%) at 2 years = cumulative number of revisions / (N - # deaths - # lost to follow up).

² Percent with Good or Excellent HSS or KSS knee score at 2 years.

³ Based on HSS knee scoring system.

⁴ Based on KSS knee scoring system.

Biomet® Oxford Partial Knee Implants in the Magnetic Resonance (MR) Environment

Biomet® Oxford Partial Knee implants are composed of non-ferromagnetic materials such as Titanium (Ti-6Al-4V), Cobalt Chrome (Co-Cr-Mo), Tantalum and ultra high molecular weight polyethylene (UHMWPE).

Biomet has performed bench testing and numerical simulations on Oxford Partial Knee implants in a Magnetic Resonance Imaging (MRI) environment. These tests determined the non-clinical effects of MRI based on scientifically relevant characteristics of the Oxford components.

MR Safety Information

The Oxford Partial Knee is determined to be MR Conditional in accordance to ASTM F2503-08 Standard Practice for Marking Devices and Other Items for Safety in the Magnetic Resonance Environment. MR Conditional refers to an item that has been demonstrated to pose no known hazards in a specified MR environment with specified conditions of use, as stated below.

Safety information for the use of MRI procedures (i.e. imaging, angiography, functional imaging, spectroscopy, etc.) pertains to shielded MRI systems under the following specifications:

- Static magnetic field of 1.5-Tesla (1.5T) and 3.0-Tesla (3.0T)
- Spatial gradient field of 2500-Gauss/cm or less
- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning
- Normal mode of the MR system
- The effects of using MR conditions above these levels have not been determined
- The effects of local RF transmit coils have not been tested and are not recommended in the area of the implant
- Per good clinical practices, the patient's legs should not touch each other and arm and legs should not touch the side of the bore during scanning

MR Information

1.5T MR system

64 MHz, GE Signa whole body coil from GE Signa MR System, model number 46-258170G1 and serial number 10146MR9; the magnet producing the static field was not present. RF power was applied continuous wave (CW) with an HP8640B (serial number 1716A06448) preamplifier and an ENI power amplifier (model 3200L, serial number 469).

- A temperature rise of 1.5 °C or less was measured at a maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of RF power application.

3.0T MR system

128MHz, GE Signa HDx 3T, Software Version=14\lx\mr Software release.14.0.M5A.0828.b, General Electric Healthcare, Milwaukee, WI. active-shielded, horizontal field scanner.

- A temperature rise of 1.4 °C or less was measured at a maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of RF power application.

Image Artifacts

MR image quality may be compromised if the area of interest is in the same area or relatively close to the position of the device.

Distortion extended as much as 8 cm from the implant in image distortion tests performed according to ASTM F2119 in a 3.0 T GE Signa HDx MR system (Software Version 14\lx\mr Software release.14.0.M5A.0828.b). Therefore, it may be necessary to optimize MR imaging parameters for the presence of these implants.

Other: Testing indicated no known risks of magnetically induced displacement force or torque.

STERILITY

Prosthetic components are sterilised by exposure to a minimum dose of 25 kGy of gamma radiation. Single Use Only. Do Not Reuse.

Do not resterilize. Do not use any component from an opened or damaged package. Do not use implants after expiration date.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

Manufactured and Distributed by:

Biomet UK Ltd.
Waterton Industrial Estate
Bridgend CF31 3XA, UK

Distributed in the United States by:

Biomet Orthopedics,
56 East Bell Drive
P.O.BOX 587
Warsaw, IN 46581 USA

CE mark on the package insert (IFU) is not valid unless there is a CE Mark on the product label.

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Biomet UK, Ltd.